

## Clinical Processing – Pictogram Key

Pictogram	Description
	Disassembling
	Assembling
	Note for assembling (written comments, reference for article-numbers is recommended)
	When a symbol is marked with a cross, the <u>action is forbidden!</u>
	Caution Magnetic! (Iron particles can build residue on the surface)
	Wrench
	Screwdriver
	Phillips screwdriver
	Stardrive® / Torx®
	Allen wrench
	Limited rectilinear motion
	Movement in one direction (Pull/Push)
	Direction of continuous rotation (Turn)
	Hold tight
	Do not bend a cable or hose
	Submerging allowed / recommended
	Ultrasonic bath allowed /recommended (Frequency 35 till 45 KHz)
	Flushing (brushing) lumina highly recommended

## Clinical Processing – Pictogram Key

	Manual Cleaning
	Automated Cleaning with manual pre-cleaning and ultrasonic cleaning
	Alcoholic wipe-disinfection (70% IPA) allowed / recommended (material compatibility according to / and disinfectant manufacturer recommendations have to be regarded)
	Alcoholic spray-disinfection (70% IPA) allowed / recommended (material compatibility according to / and disinfectant manufacturer recommendations have to be regarded)
	Steam Sterilization Detailed instructions for steam sterilization can be found in the Synthes general leaflet "Important Information"
	Do not sterilize!
	Visual inspection and Function Control Check
	Lubrication with Synthes oil (519.970 should be used)
	Lock
	Unlock
	Indicates where oil is to be applied
	Defective device. Dispose of device.
	Serviceable device